

sgRNA counting pipeline I

First things first, sharing is an excellent way to tap everyone's potential by lowering the barriers.

Authors: **Haitao Wang**, **Jimmy Zeng**, **Lakhansing Pardeshi**

Step one: Software installation

Description:

[cutadapt](#) for cutting **barcode** or **adaptor** in .fastq or .fq file

[mageck](#) for counting **activation/inhibition/knockout** sgRNA reads

```
# Before install cutadapt you should install Conda, Conda can help you install bio
software automatically, to make life easy.
# Conda installation
https://conda.io/docs/user-guide/install/macos.html#install-macos-silent

wget http://repo.continuum.io/miniconda/Miniconda3-3.7.0-Linux-x86_64.sh -O
~/miniconda.sh
bash ~/miniconda.sh -b -p $HOME/miniconda
export PATH="$HOME/miniconda/bin:$PATH"

conda config --add channels bioconda
# Cutadapt installation
conda install -c bioconda cutadapt
# MAGECK installation
conda install -c bioconda mageck

# test whether installed or not
cutadapt -h # help info to check useage
mageck -h # help info to check useage
```

Step two: Demultiplex

Description:

based on **index** or **barcodes** to deconstruct from a mix fastq file

`demultiplex.pl` this tool created by [Lakhansing Pardeshi](#) from [Chris lab](#)

Usage

```
perl demultiplex.pl --barcodes barcodes.txt --1 R1.fq --2 R2.fq --suffix <suffix to add>
```

prepare barcodes.txt look like this

```
TAAGTAGAG    HK112
ATACACGATC   HK113
GATCGCGCGGT  HK114
CGATCATGATCG    HK115
TCGATCGTTACCA   HK116
ATCGATTCCTTGGT  HK117
```

Step three trimme each fastq file for future sg counting

Description:

LentiCRISPRv2 All plasmids have the same overhangs after BsmBI digestion and the same oligos can be used for cloning into lentiCRISPRv2, lentiGuide-Puro or lentiCRISPRv1.

5' arm:GGACGAAACACCG **20bp-sg-sequence** GTTTTAGAGCTAG..... 3' arm

Target Guide Sequence Cloning Protocol

In order to clone the target sequence into the lentiCRISPRv2 or lentiGuide-Puro backbone, synthesize two oligos of the following form. All plasmids have the same overhangs after BsmBI digestion and the same oligos can be used for cloning into lentiCRISPRv2, lentiGuide-Puro or lentiCRISPRv1.

Oligo 1 → 5' — CACCGNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN — 3'
 3' — CNNNNNNNNNNNNNNNNNNNNNNNNNNNNNNCAAA — 5' ← Oligo 2

Target Sequence: 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 NGG PAM

```
# Check the CRISPR-Cas9 seq reads using 5' and 3' sequence beside sg sequence
grep --color=auto GGACGAAACACC CC1-F01_1_R1.fastq | head
grep --color=auto GTTTTAGAGCTAG CC1-F01_1_R1.fastq | head
```

Trimming seq CRISPR-Cas9 knockout screening (**usually we only use R1 in sg sequencing data**)

(compare R1 and R2, R1 have more target info than R2.)

Before Trimmed

@E00492:222:HC73FCCXY:1:1101:16133:1977 1:N:0:CTCTACTT
TC TTGTGGAAGGACGAAACACC GTCTTTGTAAGGTTCCCGTTT GTTTTAGAGCTAG AAAATAGCAAGTTAAAATAAAGGCTAGTCCGTTATCAACTTG/
+
7[F]J]]FAJ]]7J]-AAJFJ]]-AJ]]FFJ-<]]]]]]F]]]]AF-FFJ]]]]F]]]]A]]]]]]FJF]]<]]]]<A<<F-AAJF]]]]F]]AFAFJ]]]]]]AA]]]]]]F]]FF-A7-7AFA
@E00492:222:HC73FCCXY:1:1101:7202:2047 1:N:0:CTCTACTT
TC TTGTGGAAGGACGAAACACC GAGTCAAACCACTTCCCGATG GTTTTAGAGCTAG AAAATAGCAAGTTAAAATAAAGGCTAGTCCGTTATCAACTTG/
+
-]]]]]]F]]F]]FJ-<F]]]]]]]]]]7F]]7F]<AJ]-AJ]]F]]]]F]]]]]]]]]]F]]]]]]<A7F]]]]]]]]-AJ]]]]F<]]]]]]FJ]<FJF]]]]]]FJ-7FF<FFJ
@E00492:222:HC73FCCXY:1:1101:7172:2065 1:N:0:CTCTACTT
TC TTGTGGAAGGACGACACACC CGCCAAGAATTCAATTGAAA GTTTTAGAGCTAG AAAATAGCAAGTTAAAATAAAGGCTAGTCCGTTATCAACTTG/
+

After Trimmed

@E00492:222:HC73FCCXY:1:1101:16133:1977 1:N:0:CTCTACTT
GTCCTTGTAAGGTTCCCGTTT
 +
 JJ-AJJFFJF-<JJJJJJFFJ
 @E00492:222:HC73FCCXY:1:1101:7202:2047 1:N:0:CTCTACTT
GAGTCAAACCACTTCCCGATG
 +
 JJJJJJJJJ7FJJJ7FJ<AJ
 @E00492:222:HC73FCCXY:1:1101:7172:2065 1:N:0:CTCTACTT
GCGCCAAGAATTCAATTGAA
 +

```
# Trimming seq CRISPR-Cas9 knockout screening
### CRISPR-Cas9 knockout system (All plasmids have the same overhangs after BsmBI
digestion and the same oligos can be used for cloning into lentiCRISPRv2,
lentiGuide-Puro or lentiCRISPRv1)

### fq.R1 TCTTGTGGAAAGGACGAAACACCG - xxxxxxxxxxxxxxx -
GTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCTTTTTT
GAATTCGCTAGCTAG

### fq.R2
CTAGCTAGCGAATTCAAAAAGCACCGACTCGGTGCCACTTTTTCAAGTTGATAACGGACTAGCCTTATTTAACTTGCTAT
TTCTAGCTCTAAAC - xxxxxxxxxxxxxxx - CGGTGTTTCGTCCTTTCCACAAGA

# cutadapt to cut and trimme the data from 5'(only R1)
cd /PATH/
cutadapt -f fastq -q 10 \
-g GGACGAAACACC \
-o sample_1_trimmed.fastq.gz \
/PATH/sample_R1.fastq # only R1

# cutadapt to cut and trimme the data from 3' (only R1)
cutadapt -f fastq -q 10 \
-a GTTTGTAGAGCTAG \
-o /PATH/sample_2_trimmed.fastq.gz \
/PATH/sample_1_trimmed.fastq.gz

cat /PATH/sample_2_trimmed.fastq.gz > /PATH/sample_trimmed.fastq.gz
rm /PATH/sample_1_trimmed.fastq.gz
rm /PATH/sample_2_trimmed.fastq.gz
```

Step four: sg counting

Method I:

1. Raw Count sgRNA

`fastqgz_to_counts.py` this tool created by [mhorlbeck](#)

Usage (note: under python2.7; have to put required input files order like: 1. library; 2. output-path, 3. fastq files)

```
python fastqgz_to_counts.py --trim_start 1 --trim_end 21 <library.fasta> <output/path>
<all/sg.fastq.gz>
```

library

```

>0610007P14Rik_TCCTGAATGTGTTACGAAGC
tcctgaatgtgttacgaagc
>0610007P14Rik_GGTCGGGCTCCGGTACCTAG
ggtcgggctccggtacctag
>0610007P14Rik_GCCAGCTTCGTAACACATTC
gccagcttcgtaacacattc
>0610009B22Rik_TCATCATGCTGCATGACGTG
tcatcatgctgcatgacgtg
>0610009B22Rik_ATTCAATGAGTGGTTCGTCT
attcaatgagtggttcgtct
>0610009B22Rik_TCGTCACGGCTGGGCACATG
tcgtcacggctgggcacatg
>0610009D07Rik_TACACTCTGATTTGACGAAT
tacactctgatttgacgaat

```

Run

```

# count sgRNA python fastqgz_to_counts.py library-fasta output-path sg-fastq
python fastqgz_to_counts.py -p 16 \
--trim_start 1 --trim_end 21 \
GeckoV2_MGLib_A.fasta \
/output \
*trimmed.fastq.gz

```

Output

```

# count file have two columns, symbol and count
0610009020Rik_CTGTGCCAAGAGCGTTCAGC 0
0610009020Rik_TGGGTTTGGGCGTTATCCCA 0
0610010B08Rik_CGTGCATGTGAACCTCACTC 22
0610010B08Rik_GACTTCTAGAAGTTTGAAAA 6
0610010B08Rik_TATAGCTGTGAGATTCTTAT 0

```

2. Raw Count merge to table

`merged_table.pl` this tool created by [Jimmy](#) and [haitao](#)

Usage (note: can only merge files which include same row names and numbers)

```
perl merged_table.pl <file1> <file2> <file3>...
```

```
# merge selected files
perl merged_table.pl <file1> <file2> <file3>
# or All files
perl merged_table.pl /output/*
```

Output

```
# count file
sgRNA  Gene      CC1-F01 CC1-F02 CC1-F03 CC1-F04 CC1-F05
MGLibA_57638  Usp50    3    2   245  15    4
MGLibA_2481  Ablim1  10   12    5   18    6
MGLibA_18840  Foxo6   14   24   24   98   32
MGLibA_20534  Gm13040 1    505    3    7    2
MGLibA_42555  Prame    1    0    2    5    0
MGLibA_10739  Clec3b   18  621   32  122   38
MGLibA_36174  Olfr127  9    16   11   61   23
MGLibA_58314  Vmn1r238 10    6    6   69   23
MGLibA_24769  Hp1bp3   2    5    5   19    3
```

3. Normalize counts table (by total counts - size factor)

```
# add one more column named 'gene' for annotation
# sgRNA Gene sample1 sample2 ...
mageck count \
-k merged.txt \
-n /output/merged \
--norm-method total
```

Output

Normalized count file by size factor

sgRNA	Gene	CC1-F01	CC1-F02	CC1-F03	CC1-F04	CC1-F05
MGLibA_57638	Usp50	13.88454104	5.200697141	893.8208094	9.225439438	8.905259711
MGLibA_02481	Ablim1	46.28180348	31.20418285	18.24124101	11.07052733	13.35788957
MGLibA_18840	Foxo6	64.79452487	62.40836569	87.55795684	60.27287099	71.24207769
MGLibA_20534	Gm13040	4.628180348	1313.176028	10.9447446	4.305205071	4.452629855
MGLibA_42555	Prame	4.628180348	0	7.296496403	3.075146479	0
MGLibA_10739	Clec3b	83.30724627	1614.816462	116.7439425	75.0335741	84.59996725
MGLibA_36174	Olfr127	41.65362313	41.60557713	40.13073022	37.51678705	51.20524334
MGLibA_58314	Vmn1r238	46.28180348	15.60209142	21.88948921	42.43702141	51.20524334
MGLibA_24769	Hp1bp3	9.256360696	13.00174285	18.24124101	11.68555662	6.678944783

Method II:

Using `mageck one step` `count` function to count, merge and normalize data

Usage

```
mageck count -l library.txt -n <output> --sample-label name1,name2,name3 --fastq  
01_trimmed.fastq 02_trimmed.fastq 03_trimmed.fastq
```

Library

MGLibA.txt mouse library A looks like:

MGLibA_00001	TCCTGAATGTGTTACGAAGC	0610007P14Rik
MGLibA_00002	GGTCGGGCTCCGGTACCTAG	0610007P14Rik
MGLibA_00003	GCCAGCTTCGTAACACATTC	0610007P14Rik
MGLibA_00004	TCATCATGCTGCATGACGTG	0610009B22Rik
MGLibA_00005	ATTCAATGAGTGGTTCGTCT	0610009B22Rik

Run

```
mageck count -l MGLibA.txt \  
-n /Users/haitao/Desktop/sgRNA/fq \  
--sample-label CC1-F01,CC1-F02,CC1-F03 \  
--fastq CC1-F01_trimmed.fastq CC1-F02_trimmed.fastq CC1-F03_trimmed.fastq
```

Output

Normalized count file by size factor

sgRNA	Gene	CC1-F01	CC1-F02	CC1-F03	CC1-F04	CC1-F05
MGLibA_57638	Usp50	13.88454104	5.200697141	893.8208094	9.225439438	8.905259711
MGLibA_02481	Ablim1	46.28180348	31.20418285	18.24124101	11.07052733	13.35788957
MGLibA_18840	Foxo6	64.79452487	62.40836569	87.55795684	60.27287099	71.24207769
MGLibA_20534	Gm13040	4.628180348	1313.176028	10.9447446	4.305205071	4.452629855
MGLibA_42555	Prme	4.628180348	0	7.296496403	3.075146479	0
MGLibA_10739	Clec3b	83.30724627	1614.816462	116.7439425	75.0335741	84.59996725
MGLibA_36174	Olfr127	41.65362313	41.60557713	40.13073022	37.51678705	51.20524334
MGLibA_58314	Vmn1r238	46.28180348	15.60209142	21.88948921	42.43702141	51.20524334
MGLibA_24769	Hp1bp3	9.256360696	13.00174285	18.24124101	11.68555662	6.678944783

Step five: group sample comparison

Description:

Comparison between samples

MAGeCK has different commands:

`test` (if you already have count tables)

`count` (if you want to generate count tables from fastq files)

`run` (combine both test and count)

`pathway` (if you want to do the pathway test)

count file: sample.txt

sgRNA	Gene	initial1	initial2	final1	final2
A1CF_m52595977	A1CF	213	274	883	175
A1CF_m52596017	A1CF	294	412	1554	1891
AAAS_m53714382	AAAS	704	671	799	1426
AAAS_m53715169	AAAS	651	627	797	1690
AAAS_m53715176	AAAS	545	89	392	664
AAK1_m69870049	AAK1	364	465	693	2006
AATF_m35306444	AATF	449	456	1396	1402
AATF_m35306475	AATF	493	612	1102	537

Run


```

mageck test \
# Raw Count tables
-k sample.txt \
# Treatment sample labels
-t final1,final2 \
# Control sample labels
-c initial1,initial2 \
-n Output # Output labels

```

Output

id	num	neg score	neg p-value	neg fdr	neg rank	neg goodsgr	neg lfc	pos score	pos p-value	pos fdr	pos rank	pos goodsgr	pos lfc
ACRC	10	0.15729	0.38703	0.999736	31	2	0.20117	0.0026435	0.015115	0.544989	1	5	0.20117
AGAP3	10	0.95257	0.97204	0.999736	93	1	0.37513	0.0027042	0.015364	0.544989	2	8	0.37513
AGTPBP1	10	0.98731	0.98765	0.999736	98	1	0.043765	0.0049937	0.026891	0.544989	3	3	0.043765
ADCK3	10	0.99029	0.99048	0.999736	99	0	0.33152	0.0089759	0.043518	0.544989	4	7	0.33152
ABCB8	10	0.50055	0.75373	0.999736	62	3	0.3059	0.0090408	0.043797	0.544989	5	6	0.3059
ADNP	10	0.91389	0.95979	0.999736	91	2	0.35375	0.0093506	0.045005	0.544989	6	7	0.35375
ADCK1	10	0.98031	0.98274	0.999736	95	1	0.33547	0.009402	0.045324	0.544989	7	6	0.33547
ADRBK1	10	0.69259	0.86267	0.999736	76	3	0.37279	0.0095426	0.045803	0.544989	8	7	0.37279
ADCK4	10	0.57765	0.79975	0.999736	68	1	0.23437	0.013825	0.062171	0.544989	9	7	0.23437
ADRA1A	10	0.709	0.87176	0.999736	78	3	-0.11629	0.014914	0.066532	0.544989	10	3	-0.11629
ADK	10	0.5259	0.77255	0.999736	66	2	0.51714	0.016258	0.071692	0.544989	11	6	0.51714
ACTR1A	10	0.035637	0.13577	0.754297	18	5	-0.022494	0.018581	0.080384	0.544989	12	4	-0.022494
AFF4	10	0.7106	0.87272	0.999736	79	2	0.30394	0.020531	0.08733	0.544989	13	7	0.30394
ACTN4	10	0.61935	0.82197	0.999736	70	3	0.21851	0.020886	0.088568	0.544989	14	7	0.21851
AHRR	9	0.52123	0.74707	0.999736	65	3	-0.1102	0.022298	0.089306	0.544989	15	3	-0.1102
ADCK5	10	0.5362	0.77792	0.999736	67	3	0.03666	0.023019	0.096203	0.544989	16	4	0.03666
ACVR1C	10	0.99967	0.99974	0.999736	100	0	0.28521	0.023123	0.096612	0.544989	17	7	0.28521
ADARB2	10	0.91741	0.96162	0.999736	92	1	0.3507	0.024715	0.10195	0.544989	18	6	0.3507
ADRBK2	10	0.85038	0.93177	0.999736	86	2	0.28103	0.025124	0.10355	0.544989	19	7	0.28103
AAK1	10	0.23886	0.50492	0.999736	43	3	-0.025982	0.031147	0.12419	0.61874	20	3	-0.025982
AEN	10	0.039603	0.14659	0.771536	19	6	-0.31182	0.034488	0.13275	0.61874	21	1	-0.31182
AHNAK2	10	0.69378	0.86322	0.999736	77	2	0.3327	0.035849	0.13612	0.61874	22	6	0.3327
AHNAK	10	0.48698	0.74267	0.999736	61	2	0.31618	0.043452	0.15664	0.657013	23	5	0.31618
A1CF	10	0.32599	0.60137	0.999736	50	4	0.066869	0.044143	0.15841	0.657013	24	4	0.066869
AATK	10	0.6441	0.83545	0.999736	73	2	0.27164	0.047202	0.16645	0.657013	25	6	0.27164
ACVR1	10	0.98608	0.98657	0.999736	97	1	0.1924	0.049024	0.17082	0.657013	26	5	0.1924
ACSS2	10	0.88811	0.94737	0.999736	89	1	0.18848	0.051892	0.17813	0.659736	27	6	0.18848

Step Six*: sg downstream analysis

coming soon